

Rab5 Rabbit mAb


ABZ5032

Product Overview


Name	Rab5 Rabbit mAb
Accession(Primary)	P20339
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 24 kDa; Observed MW: 24 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200

Target information (P20339)

Gene Name	RAB5A
About Gene Symbol	RAB5A (RAB5A, Member RAS Oncogene Family) is a small GTPase involved in cell signaling and proliferation. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Rab5 in MCF-7 lysates using Rab5 antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon using Rab5 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.