

Hsp90 alpha/beta Rabbit mAb


ABZ5034

Product Overview


Name	Hsp90 alpha/beta Rabbit mAb
Accession(Primary)	P07900
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 85 kDa; Observed MW: 90 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:1000; ICC/IF-1:100-1:200; FC-1:200-1:500; IP-1:50-1:100

Target information (P07900)


Gene Name	HSP90AA1/HSP90AB1
About Gene Symbol	



Western blot analysis of Hsp90 in (1) HeLa lysates; (2) Jurkat lysates; (3) RAW264.7 lysates; (4) NIH/3T3 lysates; (5) PC-12 lysates; (6) C6 lysates using Hsp90 alpha/beta antibody.



Immunohistochemistry analysis of paraffin-embedded Human liver carcinoma using Hsp90 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunofluorescence analysis of Hsp90 alpha/beta in HeLa using Hsp90 antibody.