

# Histone H3.3 Rabbit mAb


ABZ5035

## Product Overview


|                             |                                                             |
|-----------------------------|-------------------------------------------------------------|
| <b>Name</b>                 | Histone H3.3 Rabbit mAb                                     |
| <b>Accession(Primary)</b>   | P84243                                                      |
| <b>Host Species</b>         | Rabbit                                                      |
| <b>Clonality</b>            | Monoclonal                                                  |
| <b>Isotype</b>              | IgG                                                         |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                                           |
| <b>Calculated MW</b>        | Calculated MW: 15 kDa; Observed MW: 15 kDa                  |
| <b>Recommended Dilution</b> | WB-1:500-1:1000; IHC-P-1:200-1:1000;<br>ICC/IF-1:500-1:1000 |

## Target information (P84243)

|                          |       |
|--------------------------|-------|
| <b>Gene Name</b>         | H3-3A |
| <b>About Gene Symbol</b> |       |



Western blot analysis of Histone H3.3 in (1) HeLa lysates; (2) NIH/3T3 lysates using Histone H3.3 antibody.



Immunohistochemistry analysis of paraffin-embedded mouse stomach using Histone H3.3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.