

Midkine Rabbit mAb


ABZ5036

Product Overview


Name	Midkine Rabbit mAb
Accession(Primary)	P21741
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 16 kDa; Observed MW: 16 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:20-1:50; ICC/IF-1:50-1:100; FC-1:100-1:200; IP-1:20

Target information (P21741)

Gene Name	MDK
About Gene Symbol	



Western blot analysis of Midkine in Midkine Recombinant protein using Midkine antibody.



Immunohistochemistry analysis of paraffin-embedded Human liver carcinoma using Midkine antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.