

Prion Protein PrP Rabbit mAb


ABZ5042

Product Overview

Name	Prion Protein PrP Rabbit mAb
Accession(Primary)	P04156
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 28 kDa; Observed MW: 28 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500; FC-1:20-1:50

Target information (P04156)

Gene Name	PRNP
About Gene Symbol	



Western blot analysis of Prion Protein in Hunam fetal brain lysates using Prion Protein antibody.