


# Cytochrome C Rabbit mAb

ABZ5043


## Product Overview

|                             |                                                 |
|-----------------------------|-------------------------------------------------|
| <b>Name</b>                 | Cytochrome C Rabbit mAb                         |
| <b>Accession(Primary)</b>   | P99999                                          |
| <b>Host Species</b>         | Rabbit                                          |
| <b>Clonality</b>            | Monoclonal                                      |
| <b>Isotype</b>              | IgG                                             |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                               |
| <b>Calculated MW</b>        | Calculated MW: 12 kDa; Observed MW: 15 kDa      |
| <b>Recommended Dilution</b> | WB-1:500-1:1000; IHC-1:100-1:200; IF-1:50-1:200 |


## Target information (P99999)

**Gene Name** CYCS**About Gene Symbol**

Western blot analysis of Cytochrome C in HeLa lysates using Cytochrome C antibody.



Immunohistochemistry analysis of paraffin-embedded Human colonic cancer using Cytochrome C antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunofluorescence analysis of Cytochrome C in Hela using Cytochrome C antibody.