

Radixin Rabbit mAb


ABZ5044

Product Overview


Name	Radixin Rabbit mAb
Accession(Primary)	P35241
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 69 kDa; Observed MW: 80 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:20-1:50; FC-1:50-1:100

Target information (P35241)

Gene Name	RDX
About Gene Symbol	RDX (Radixin) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in hearing loss and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Radixin in HeLa lysates using Radixin antibody.



Immunohistochemistry analysis of paraffin-embedded Human spleen using Radixin antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.