

Rb Rabbit mAb


ABZ5045

Product Overview


Name	Rb Rabbit mAb
Accession(Primary)	P06400
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 106 kDa; Observed MW: 106 kDa
Recommended Dilution	WB-1:500-1:1000; IHC:1:200-1:1000; ICC/IF-1:200-1:500; FC-1:50-1:100; IP-1:20-1:50

Target information (P06400)


Gene Name	RB1
About Gene Symbol	RB1 (RB Transcriptional Corepressor 1) is a protein phosphatase involved in gene transcription and regulation. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Retinoblastoma in (1) Jurkat lysates; (2) MCF-7 lysates using Rb antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using Retinoblastoma antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunofluorescence analysis of Rb in SH-SY5Y using Retinoblastoma antibody.