

Asparagine Synthetase Rabbit mAb

ABZ5049


Product Overview

Name	Asparagine Synthetase Rabbit mAb
Accession(Primary)	P08243
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 64 kDa; Observed MW: 61 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100


Target information (P08243)

Gene Name ASNS

About Gene Symbol



Western blot analysis of ASNS in K562 lysates using Asparagine Synthetase antibody.



Immunohistochemistry analysis of paraffin-embedded Human bladder using ASNS antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.