

MEKK2 Rabbit mAb


ABZ5050

Product Overview


Name	MEKK2 Rabbit mAb
Accession(Primary)	Q9Y2U5
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 70 kDa; Observed MW: 70 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; ICC/IF-1:200-1:500; FC-1:50-1:100; IP-1:20-1:50

Target information (Q9Y2U5)

Gene Name	MAP3K2
About Gene Symbol	MAP3K2 (Mitogen-Activated Protein Kinase Kinase 2) is a mitogen-activated protein kinase involved in MAPK/ERK signaling. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of MEKK2 in HeLa lysates using MEKK2 antibody.



Immunohistochemistry analysis of paraffin-embedded Human thyroid cancer using MEKK2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.