

VAV2 Rabbit mAb


ABZ5051

Product Overview


Name	VAV2 Rabbit mAb
Accession(Primary)	P52735
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 101 kDa; Observed MW: 101 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; FC-1:200-1:1000; IP-1:20

Target information (P52735)

Gene Name	VAV2
About Gene Symbol	



Western blot analysis of VAV2 in 293T lysates using VAV2 antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using VAV2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.