

Topoisomerase I Rabbit mAb


ABZ5053

Product Overview


Name	Topoisomerase I Rabbit mAb
Accession(Primary)	P11387
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 91 kDa; Observed MW: 91 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:200-1:500; FC-1:10-1:100

Target information (P11387)

Gene Name	TOP1
About Gene Symbol	



Western blot analysis of TOP1 in Jurkat lysates using Topoisomerase I antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon using TOP1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.