

INPP4B Rabbit mAb


ABZ5055

Product Overview

Name	INPP4B Rabbit mAb
Accession(Primary)	O15327
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 105 kDa; Observed MW: 105 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:20-1:50; ICC/IF-1:20-1:50; FC-1:50-1:100

Target information (O15327)

Gene Name	INPP4B
About Gene Symbol	



Western blot analysis of INPP4B in MCF-7 lysates
using INPP4B antibody.