

IL-1RA Protein Rabbit mAb

ABZ5057


Product Overview

Name	IL-1RA Protein Rabbit mAb
Accession(Primary)	P18510
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 20 kDa; Observed MW: 20 kDa
Recommended Dilution	WB-1:2000-1:5000

Target information (P18510)

Gene Name IL1RN

About Gene Symbol



Western blot analysis of IL-1RA in (1) HeLa lysates; (2) 3T3 lysates using IL-1 Receptor Antagonist Protein antibody.

