


SMC1A Rabbit mAb

ABZ5059

Product Overview

Name	SMC1A Rabbit mAb
Accession(Primary)	Q14683
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 143 kDa; Observed MW: 143 kDa
Recommended Dilution	WB-1:2000-1:5000; IHC-P-1:200-1:500

Target information (Q14683)

Gene Name SMC1A**About Gene Symbol**

Western blot analysis of SMC1 in (1) HeLa lysates;
(2) Jurkat lysates using SMC1A antibody.

