

Histone H2B Rabbit mAb


ABZ5060

Product Overview


Name	Histone H2B Rabbit mAb
Accession(Primary)	O60814
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 14 kDa; Observed MW: 14 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC:1:20; ICC/IF-1:20; FC-1:20; IP-1:20

Target information (O60814)

Gene Name	H2BC12
About Gene Symbol	



Western blot analysis of Histone H2B in A431 lysates using Histone H2B antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using Histone H2B antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.