

MEK3 Rabbit mAb


ABZ5062

Product Overview


Name	MEK3 Rabbit mAb
Accession(Primary)	P46734
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 39 kDa; Observed MW: 39 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500; ICC/IF-1:100-1:200; FC-1:200-1:1000; IP-1:50-1:100

Target information (P46734)

Gene Name	MAP2K3
About Gene Symbol	MAP2K3 (Mitogen-Activated Protein Kinase Kinase 3) is a mitogen-activated protein kinase involved in MAPK/ERK signaling. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of MEK3 in HepG2 lysates using MEK3 antibody.



Immunohistochemistry analysis of paraffin-embedded Human liver using MEK3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.