

HDAC8 Rabbit mAb


ABZ5065

Product Overview

Name	HDAC8 Rabbit mAb
Accession(Primary)	Q9BY41
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 42 kDa; Observed MW: 42 kDa
Recommended Dilution	WB-1:1000-1:5000; FC-1:50-1:100; IP-1:20-1:50

Target information (Q9BY41)

Gene Name	HDAC8
About Gene Symbol	HDAC8 (Histone Deacetylase 8) is a histone deacetylase involved in chromatin remodeling and epigenetics. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of HDAC8 in HeLa lysates
using HDAC8 antibody.