

TRAP1 Rabbit mAb


ABZ5069

Product Overview

Name	TRAP1 Rabbit mAb
Accession(Primary)	Q12931
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 80 kDa; Observed MW: 75 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (Q12931)

Gene Name	TRAP1
About Gene Symbol	TRAP1 (TNF Receptor Associated Protein 1) is a molecular chaperone involved in immune signaling. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Hsp75 in K562 lysates
using Hsp75 antibody.