


HspA2 Rabbit mAb

ABZ5071

Product Overview

Name	HspA2 Rabbit mAb
Accession(Primary)	P54652
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 70 kDa; Observed MW: 70 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:2000-1:5000

Target information (P54652)

Gene Name HSPA2**About Gene Symbol**

Western blot analysis of HSPA2 in MCF-7 lysates using HspA2 antibody.

