

JNK1/3 Rabbit mAb


ABZ5072

Product Overview

Name	JNK1/3 Rabbit mAb
Accession(Primary)	P45983
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 48 kDa; Observed MW: 46,54 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20

Target information (P45983)

Gene Name	MAPK8/MAPK10
About Gene Symbol	MAPK8 (Mitogen-Activated Protein Kinase 8) is a mitogen-activated protein kinase involved in MAPK/ERK signaling. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of JNK1/JNK3 in HeLa lysates
using JNK1/3 antibody.