

KMT6/EZH2 Rabbit mAb


ABZ5075

Product Overview


Name	KMT6/EZH2 Rabbit mAb
Accession(Primary)	Q15910
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 85 kDa; Observed MW: 98 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IF-1:50-1:200; FC-1:50-1:100

Target information (Q15910)


Gene Name	EZH2
About Gene Symbol	EZH2 (Enhancer Of Zeste 2 Polycomb Repressive Complex 2 Subunit 2) is a methyltransferase involved in chromatin remodeling and epigenetics. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of KMT6 / EZH2 in HEK293 lysates using KMT6 antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using KMT6 / EZH2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunofluorescence analysis of KMT6 in HeLa using KMT6 / EZH2 antibody.