

Glucocorticoid Receptor Rabbit mAb


ABZ5081

Product Overview

Name	Glucocorticoid Receptor Rabbit mAb
Accession(Primary)	P04150
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 86 kDa; Observed MW: 95 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:2000; ICC/IF-1:200-1:500; FC-1:200-1:500

Target information (P04150)

Gene Name	NR3C1
About Gene Symbol	NR3C1 (Nuclear Receptor Subfamily 3 Group C Member 1) is a membrane receptor involved in receptor-mediated signaling. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mec



Western blot analysis of GR in A431 lysates using
Glucocorticoid Receptor antibody.