

BMP11 Rabbit mAb


ABZ5083

Product Overview

Name	BMP11 Rabbit mAb
Accession(Primary)	O95390
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 45 kDa; Observed MW: 45 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:500-1:2000; ICC/IF-1:100-1:200; FC-1:100-1:200

Target information (O95390)

Gene Name	GDF11
About Gene Symbol	



Western blot analysis of BMP11 in (1) HeLa lysates; (2) C6 lysates using BMP11 antibody.