

Thymidine Kinase 1/TK1 Rabbit mAb


ABZ5084

Product Overview

Name	Thymidine Kinase 1/TK1 Rabbit mAb
Accession(Primary)	P04183
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 25 kDa; Observed MW: 25 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:100-1:200; FC-1:100-1:200; IP-1:20

Target information (P04183)

Gene Name	TK1
About Gene Symbol	TK1 (Thymidine Kinase 1) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Thymidine Kinase 1 in Molt4 lysates using Thymidine Kinase antibody.