


PDIA2/PDIP Rabbit mAb

ABZ5085

Product Overview

Name	PDIA2/PDIP Rabbit mAb
Accession(Primary)	Q13087
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 58 kDa; Observed MW: 58 kDa
Recommended Dilution	WB-1:1000-1:5000; FC-1:50-1:100

Target information (Q13087)

Gene Name PDIA2**About Gene Symbol**

Western blot analysis of PDIA2 in HeLa lysates using PDIA2 antibody.

