


Sarcomeric alpha Actinin Rabbit mAb

ABZ5086

Product Overview

Name	Sarcomeric alpha Actinin Rabbit mAb
Accession(Primary)	P35609
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 104 kDa; Observed MW: 104 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IP-1:20

Target information (P35609)

Gene Name ACTN2**About Gene Symbol**

Western blot analysis of ACTN2 in HepG2 lysates using Sarcomeric alpha Actinin antibody.

