

MiTF Rabbit mAb


ABZ5089

Product Overview

Name	MiTF Rabbit mAb
Accession(Primary)	O75030
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 59 kDa; Observed MW: 52 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (O75030)

Gene Name	MITF
About Gene Symbol	MITF (Melanocyte Inducing Transcription Factor) is a transcription factor involved in gene transcription and regulation. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of MiTF in A375 lysates
using MiTF antibody.