


PROX1 Rabbit mAb

ABZ5091

Product Overview

Name	PROX1 Rabbit mAb
Accession(Primary)	Q92786
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 83 kDa; Observed MW: 95 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (Q92786)

Gene Name PROX1**About Gene Symbol**

Western blot analysis of PROX1 in HepG2 lysates using PROX1 antibody.

