

MHC Class I Rabbit mAb


ABZ5092

Product Overview

Name	MHC Class I Rabbit mAb
Accession(Primary)	P30443
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 41 kDa; Observed MW: 41 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:500-1:5000; ICC/IF-1:200-1:500; FC-1:50-1:100

Target information (P30443)

Gene Name	HLA-A
About Gene Symbol	



Western blot analysis of MHC class I in Raji
lysates using MHC Class I antibody.