

Human IgA Rabbit mAb


ABZ5093

Product Overview

Name	Human IgA Rabbit mAb
Accession(Primary)	P01876
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 60 kDa; Observed MW: 60 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IF-1:50-1:200; IP-1:20

Target information (P01876)

Gene Name	IGHA1/IGHA2
About Gene Symbol	



Western blot analysis of IgA in Human plasma
lysates using Human IgA antibody.