

Calbindin Rabbit mAb


ABZ5094

Product Overview

Name	Calbindin Rabbit mAb
Accession(Primary)	P05937
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 30 kDa; Observed MW: 30 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:1000; ICC/IF-1:50-1:100; FC-1:50-1:100

Target information (P05937)

Gene Name	CALB1
About Gene Symbol	



Western blot analysis of Calbindin in Human fetal brain lysates using Calbindin antibody.