

Human kappa Light Chain Rabbit mAb


ABZ5096

Product Overview

Name	Human kappa Light Chain Rabbit mAb
Accession(Primary)	P01834
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 12 kDa; Observed MW: 25 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:1000-1:5000; IP-1:10-1:100

Target information (P01834)

Gene Name	IGKC
About Gene Symbol	



Western blot analysis of Human Kappa light chain in Human plasma lysates using Human kappa Light Chain antibody.