

GATA1 Rabbit mAb


ABZ5108

Product Overview


Name	GATA1 Rabbit mAb
Accession(Primary)	P15976
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 43 kDa; Observed MW: 43 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IF-1:50-1:200; IP-1:20; ChIP-1:20

Target information (P15976)

Gene Name	GATA1
About Gene Symbol	



Western blot analysis of GATA1 in K562 lysates using GATA1 antibody.



Immunohistochemistry analysis of paraffin-embedded mouse spleen using GATA1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.