

Profilin 1 Rabbit mAb


ABZ5109

Product Overview

Name	Profilin 1 Rabbit mAb
Accession(Primary)	P07737
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 15 kDa; Observed MW: 15 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:1000; ICC/IF-1:50-1:100; FC-1:100-1:200

Target information (P07737)

Gene Name	PFN1
About Gene Symbol	



Western blot analysis of Profilin1 in Jurkat lysates
using Profilin 1 antibody.