

PCSK9 Rabbit mAb

ABZ5112

Product Overview

Name	PCSK9 Rabbit mAb
Accession(Primary)	Q8NBP7
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 74 kDa; Observed MW: 62,78 kDa
Recommended Dilution	WB-1:1000-1:5000; IP-1:20-1:50

Target information (Q8NBP7)

Gene Name PCSK9**About Gene Symbol**

