


Glucose Transporter GLUT3 Rabbit mAb

ABZ5115

Product Overview

Name	Glucose Transporter GLUT3 Rabbit mAb
Accession(Primary)	P11169
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 54 kDa; Observed MW: 54 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (P11169)

Gene Name SLC2A3**About Gene Symbol**

Western blot analysis of GLUT3 in (1) HepG2 lysates; (2) NIH/3T3 lysates using Glucose Transporter GLUT3 antibody.

