

SNAP25 Rabbit mAb


ABZ5116

Product Overview

Name	SNAP25 Rabbit mAb
Accession(Primary)	P60880
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 23 kDa; Observed MW: 23 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200; IP-1:20

Target information (P60880)

Gene Name	SNAP25
About Gene Symbol	SNAP25 (Synaptosome Associated Protein 25) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of SNAP25 in SH-SY5Y
lysates using SNAP25 antibody.