

Human IgG Rabbit mAb


ABZ5117

Product Overview


Name	Human IgG Rabbit mAb
Accession(Primary)	P01857
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 55,25 kDa; Observed MW: 55,25 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:1000

Target information (P01857)

Gene Name	IGHG1
About Gene Symbol	



Western blot analysis of Human IgG in Human plasma lysates using Human IgG antibody.



Immunohistochemistry analysis of paraffin-embedded Human spleen using Human IgG antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.