


cIAP1 Rabbit mAb

ABZ5118

Product Overview

Name	cIAP1 Rabbit mAb
Accession(Primary)	Q13490
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 70 kDa; Observed MW: 70 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500

Target information (Q13490)

Gene Name BIRC2**About Gene Symbol**

Western blot analysis of BIRC2 in HepG2 lysates using cIAP1 antibody.

