


PLGF Rabbit mAb

ABZ5119

Product Overview

Name	PLGF Rabbit mAb
Accession(Primary)	P49763
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 25 kDa; Observed MW: 25 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (P49763)

Gene Name PGF**About Gene Symbol**

Western blot analysis of PLGF in Human PLGF
Recombinant protein using PLGF antibody.

