

Parathyroid Hormone Rabbit mAb

ABZ5120


Product Overview

Name	Parathyroid Hormone Rabbit mAb
Accession(Primary)	P01270
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 13 kDa; Observed MW: 36 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200

Target information (P01270)

Gene Name

PTH

About Gene Symbol

Western blot analysis of Parathyroid Hormone in Human Parathyroid Hormone recombinant protein using Parathyroid Hormone antibody.

