

Lambda Light Chain Rabbit mAb


ABZ5121

Product Overview

Name	Lambda Light Chain Rabbit mAb
Accession(Primary)	P01701
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 25 kDa; Observed MW: 25 kDa
Recommended Dilution	WB-1:500-1:1000; IHC:1:100-1:500; ICC/IF-1:100-1:200

Target information (P01701)

Gene Name	IGLV1-51/IGLC1
About Gene Symbol	



Western blot analysis of Lambda Light chain in Human plasma lysates using Lambda Light Chain antibody.