

Calretinin Rabbit mAb


ABZ5122

Product Overview

Name	Calretinin Rabbit mAb
Accession(Primary)	P22676
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 32 kDa; Observed MW: 29 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:500-1:2000; ICC/IF-1:100-1:200; FC-1:20-1:50

Target information (P22676)

Gene Name	CALB2
About Gene Symbol	



Western blot analysis of Calretinin in SH-SY5Y
lysates using Calretinin antibody.