

N WASP Rabbit mAb


ABZ5124

Product Overview

Name	N WASP Rabbit mAb
Accession(Primary)	O00401
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 55 kDa; Observed MW: 55 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; ICC/IF-1:100-1:200; FC-1:100-1:200

Target information (O00401)

Gene Name	WASL
About Gene Symbol	



Western blot analysis of WASL in HeLa lysates
using N WASP antibody.