

TCTP Rabbit mAb


ABZ5132

Product Overview

Name	TCTP Rabbit mAb
Accession(Primary)	P13693
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 20 kDa; Observed MW: 20 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:200-1:500; FC-1:100-1:1000

Target information (P13693)

Gene Name	TPT1
About Gene Symbol	



Western blot analysis of TCTP Receptor in HeLa lysates using TCTP antibody.