

Growth Hormone Receptor Rabbit mAb


ABZ5133

Product Overview

Name	Growth Hormone Receptor Rabbit mAb
Accession(Primary)	P10912
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 72 kDa; Observed MW: 140 kDa
Recommended Dilution	WB-1:1000-1:15000; ICC/IF-1:100-1:200; FC-1:200-1:500

Target information (P10912)

Gene Name	GHR
About Gene Symbol	



Western blot analysis of Growth hormone receptor in (1) MCF-7 lysates; (2) mouse brain lysates using Growth Hormone Receptor antibody.