

# Adenosine A1 Receptor Rabbit mAb


ABZ5134

## Product Overview

|                             |                                                                                            |
|-----------------------------|--------------------------------------------------------------------------------------------|
| <b>Name</b>                 | Adenosine A1 Receptor Rabbit mAb                                                           |
| <b>Accession(Primary)</b>   | P30542                                                                                     |
| <b>Host Species</b>         | Rabbit                                                                                     |
| <b>Clonality</b>            | Monoclonal                                                                                 |
| <b>Isotype</b>              | IgG                                                                                        |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                                                                          |
| <b>Calculated MW</b>        | Calculated MW: 37 kDa; Observed MW: 37 kDa                                                 |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000; IHC-P-1:100-1:1000;<br>ICC/IF-1:50-1:100; FC-1:100-1:1000; IP-1:20-1:100 |

## Target information (P30542)

|                          |        |
|--------------------------|--------|
| <b>Gene Name</b>         | ADORA1 |
| <b>About Gene Symbol</b> |        |



Western blot analysis of ADORA1 in A549 lysates  
using Adenosine A1 Receptor antibody.