

CD55 Rabbit mAb


ABZ5135

Product Overview

Name	CD55 Rabbit mAb
Accession(Primary)	P08174
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 41 kDa; Observed MW: 55-80 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100

Target information (P08174)

Gene Name	CD55
About Gene Symbol	CD55 (CD55 Molecule (Cromer Blood Group)) is a cell surface marker involved in adaptive immunity. It is an important research target in immunology, flow cytometry, and cancer immunotherapy and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of CD55 in A549 lysates
using CD55 antibody.