


Follistatin Rabbit mAb

ABZ5137

Product Overview

Name	Follistatin Rabbit mAb
Accession(Primary)	P19883
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Rat
Calculated MW	Calculated MW: 38 kDa; Observed MW: 38 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (P19883)

Gene Name FST**About Gene Symbol**

Western blot analysis of Follistatin in HepG2
lysates using Follistatin antibody.

